

# “Commissioning and Quality Assurance of Volumetric Modulated Arc Radiotherapy with RapidArc”

Adriana Martins, André Pereira, Mara Barreiros, Sofia Nunes

Hospital do Espírito Santo de Évora, Portugal

European Medical Physics and Engineering Conference 2011

Trinity College Dublin, Ireland



Adriana Martins

- September 2011 -

# Overview

---

1. Introduction
2. Commissioning and QA  
Tests for VMAT
3. Concluding Remarks



# 1. Introduction



## INTENSITY MODULATED RADIATION THERAPY (IMRT)

### Step & Shoot

- Static Gantry
- Static MLC
- Multiple Fields



### IMRT delivery



### Sliding Window

- Static Gantry
- Dynamic MLC
- Multiple Fields

### Tomotherapy

- Rotational fan beams
- Slice-by-slice

### IMAT

- Dynamic Gantry
- Dynamic MLC
- Multiple Fields

### VMAT

- Dynamic Gantry
- Dynamic MLC
- Variable Dose Rate
- Variable Gantry Speed
- Single Field



Reduced treatment time (2min)

# 1. Introduction

## MAJOR QA CONCERNS IN IMRT - MLC

### 3DCRT

MLC defines dose near the the borders of the critical structures.

- Affects high dose gradient regions
- Variations of leaf speed have no effect on the dose distribution
- No use of closed leaf pairs

### DYNAMIC MLC IMRT

Leaves modulate the dose delivered to target volume.

- Dose delivered is sensitive to the **gap width** defined by each pair  
Average gap error  $\leq 0,2$  mm for an average dose error  $< 1\%$   
(LoSasso et al., 03)
- **Leaf speed** may influence the dose distribution
- Effect of rounded leaf edge **transmission from closed leaves**


# 1. Introduction

## VMAT WITH RAPIDARC (Varian Medical Systems)

Both the treatment planning and the linac system incorporate in a single field/arc:

- Dynamic MLC
- Variable dose rate
- Variable gantry speed

- Optimization done using progressive sampling
- Constraints imposed on MLC motion and MU/degree



# 1. Introduction

---

## PURPOSE

- Clinical implementation of VMAT (from 3DCRT)

# 2. Commissioning and QA Tests

---

2.1. Equipment

2.2. Dynamic MLC

2.3. Dose Rate

2.4. Gantry

2.5. RapidArc QA Combined Tests (MLC, Gantry, Dose Rate)

2.6. Treatment Planning System

## 2.1. EQUIPMENT

---

### LINAC

- LINAC: Clinac DHX (Varian Medical Systems)
- Energy: 6MV

### TPS

- TPS: Eclipse (Varian Medical Systems)
- Calculation Algorithm: AAA
- Optimization Algorithm: PRO



## 2.1. EQUIPMENT

---

### DOSIMETRIC EQUIPMENT

- Ionization chamber: cylindrical 0,125cc (PTW)
- 2D array: seven29 (PTW)
- Film Dosimetry: EBT2 (Standard Imaging) with Epson Scanner 10000XL and RIT 113

### PHANTOMS

- CIRS
- Octavius (PTW)
- RW3 Slab phantom (PTW)

## 2.2. DYNAMIC MLC TESTS

---

### Mechanical

- Leaf positioning/gap width
- Leaf Speed

### Dosimetric

- Gap stability
- Dosimetric leaf gap



## 2.2. DYNAMIC MLC TESTS (LoSasso *et al*, 98, 01, 03)

---

### LEAF POSITIONING/GAP WIDTH


#### FENCE TESTS

**Setup:** Film and EPID at isocenter

- Picket Fence, Garden Fence
- Match lines of 1 mm gap
- $G=0^{\circ}/90^{\circ}/270^{\circ}$
- Introduction of different positioning errors.

**Analysis:** Visual (0.5 mm)

RIT Software (Position, FWHM , Transmission)



## 2.2. DYNAMIC MLC TESTS (LoSasso *et al*, 98, 01, 03)

### LEAF POSITIONING/GAP WIDTH

#### SWEEPING GAP:

**Setup:** 2D array in Octavius at isocenter

- 5 mm sliding window beam creating a uniform dynamic field of 12x27 cm
- Dose normalization to static field
- Center and off-axis (2cm) measurements
- $G=0^{\circ}/90^{\circ}/270^{\circ}$

**Analysis:** 2% tolerance (Van Esch *et al.*, 02)

RESULTS (comparing to $G=0^{\circ}$ )			
Average	X=0.0 cm	X=2.0 cm	X=-2.0 cm
G=90°	0,2%	0,2%	-0,1%
G=270°	1,2%	1,1%	-0,8%

## 2.2. DYNAMIC MLC TESTS

### LEAF SPEED



**Setup:** EPID at isocenter

- dMLC plan (sweeping gap 1cm) with all the leaves moving at constant speed.
- Three plans: low (0,4cm/s), medium (1,2cm/s) and high (2.4cm/s) speed.
- Field dimensions of 10x40 cm ( $C=90^\circ$ ) and 10x30 cm ( $C=0^\circ$ ).
- $G=0^\circ/90^\circ/270^\circ$

**Analysis:** Visual

Profile Uniformity

Dynalog File Viewer



## 2.2. DYNAMIC MLC TESTS (LoSasso *et al.*, 98, 01, 03)

### DOSIMETRIC LEAF GAP

**Setup:** Ionization chamber at isocenter with solid water phantom

- 6 different dMLC fields with fixed leaf gap widths (0.5; 1.0; 2.0; 5.0; 10.0; 20.0 mm) sweeping at constant velocity
- Total field size of 10x10 cm.


**Analysis:** Extrapolation to zero dose gives the DLG

#### Literature:

- 1.9 mm (Sivakumar, S.S. *et al.*, 08)
- 1.8 mm (LoSasso, 03)
- 1.5 mm (Bhardwaj *et al.*, 07)

#### RESULTS

DLG	1.7 mm
-----	--------



## 2.3. DOSE RATE TESTS

---

- Linearity
- Stability
- Symmetry and Homogeneity






## 2.3. DOSE RATE TESTS

### LINEARITY

**Setup:** Ionization chamber at isocenter in solid water phantom

- Irradiation of fields with the same size and increasing MU
- Dose Normalization to Reference Field of 100 MU
- 100 MU/min, 300MU/min and 600 MU/min.

#### **Analysis:**





## 2.3. DOSE RATE TESTS (Yang *et al*, 04)

---

### STABILITY

**Setup:** Ionization chamber at isocenter in solid water phantom

- Irradiation of different combination of (number of fields, MU) with the same total dose - from 1 field with 100 MU up to 100 fields with 1 MU.
- 100 UM/min, 300UM/min and 600 MU/min.
- Dose normalization to Reference Field of 100 MU
- Dose normalization to Reference Dose Rate of 300 MU/min

**Analysis:** 2% Tolerance

## 2.3. DOSE RATE TESTS (Yang et al, 04)

### STABILITY

RESULTS Dose Normalization to Ref Field				
UM/FIELD	# FIELDS	100 MU/min	300 MU/min	600 MU/min
100	1	0.0 %	0.0 %	0.0 %
20	5	-0.1 %	-0.1 %	-0.2 %
10	10	-0.2 %	-0.3 %	-0.5 %
5	20	-0.4 %	-0.6 %	-1.1 %
4	25	-0.4 %	-0.7 %	-1.3 %
2	50	-1.1 %	-1.4 %	-2.2 %
1	100	-2.2%	-2.8 %	-5.4 %

RESULTS Dose Normalization to Ref Dose Rate			
UM/FIELD	# FIELDS	100 MU/min	600 MU/min
100	1	0.0 %	-0.1 %
20	5	0.1 %	-0.2 %
10	10	0.1 %	-0.3 %
5	20	0.2 %	-0.6 %
4	25	0.3 %	-0.7 %
2	50	0.3 %	-1.0 %
1	100	0.7%	-2.6 %

## 2.3. DOSE RATE TESTS

### SIMMETRY AND HOMOGENEITY


**Setup:** 2D array with Octavius phantom

- Irradiation of fields of different sizes (5x5 cm, 10x10cm, 20x20 cm) and increasing MU.
- 100 MU/min, 300MU/min and 600 MU/min.

**Analysis:** 3% Tolerance

#### RESULTS

Variation of  $\pm 0.5\%$   
relative to 50 MU Field



## 2.4. GANTRY TESTS

---

### DOSE STABILITY

**Setup:** Ionization chamber at isocenter with buildup cap

- Dose comparison between  $G=0^\circ$  and  $360^\circ$  Arc rotation
- 9 Field dimensions: from 5x5 cm up to 40x40 cm.

**Analysis:** 1% Tolerance

RESULTS
$\leq 0.5\%$

## 2.5. RAPIDARC QA COMBINED TESTS

---

- Picket Fence
- Gantry speed vs Dose Rate vs Output
- MLC Speed vs Output



## 2.5. RAPIDARC QA COMBINED TESTS

(Ling *et al*, 08)

### GANTRY SPEED VS DOSE RATE VS OUTPUT

**Setup:** EPID at isocenter

- 7 static MLC fields of 20x1.8 cm irradiated with the same dose but different combinations of dose rate(111-600 MU/min) and gantry speed (5.5-4.3<sup>0</sup>/s).
- Normalization to a static open field.

#### **Analysis:**

RESULTS (Central Profile)								
Band	B1	B2	B3	B4	B5	B6	B7	MEAN
Deviation	1.9 %	0.7 %	0.6 %	-0.4%	-1.0 %	-0.8 %	-1.0 %	0.9 %

#### **Literature:**

- Mean 0.7 % (Ling *et al.*, 08)

## 2.5. RAPIDARC QA COMBINED TESTS

(Ling *et al*, 08)

### MLC SPEED VS OUTPUT

**Setup:** EPID at isocenter

- 4 bands of dynamic MLC fields combined different leaf speeds (0.46-2.76 cm/s) with different dose rates to achieve the same dose.
- Normalization to a static open field.

### **Analysis:**

RESULTS (Central Profile)					
Band	B1	B2	B3	B4	MEAN
Deviation	0.1 %	-0.2 %	0.3 %	-0.2%	0.2 %

### **Literature:**

- Mean 0.4 % (Ling *et al.*, 08)

## 2.6. TPS TESTS

---

- MLC parameter optimization
- QA RapidArc Test Plans
- Clinical RapidArc Test Plans






## 2.6. TPS TESTS (Van Esch *et al*, 02)

### MLC PARAMETER OPTIMIZATION

**Setup:** Ionization chamber in water phantom

- Chair-like fluence was created and irradiated
- Measurements at depths: 5 cm, 10 cm and 20 cm
- Measurements at: Homogeneous dose, Transmission and DLG zones
- Comparison of calculated versus measured dose for optimization of MLC transmission and DLG parameters in TPS.



## 2.6. TPS TESTS (Van Esch *et al*, 02)

### MLC PARAMENTER OPTIMIZATION

**Analysis:** 3% Tolerance

RESULTS			
DEPTH	ZONE	INITIAL	FINAL
5 cm	ISO	-0.8 %	-0.4 %
	TRANS	-10.2 %	0.8 %
	DLG	-7.6 %	-0.6 %
10 cm	ISO	-0.6 %	-0.3 %
	TRANS	-9.7 %	-0.7 %
	DLG	-5.3 %	-0.5 %
20 cm	ISO	-0.8 %	-0.6 %
	TRANS	-2.8 %	-1.8 %
	DLG	-1.2 %	1.6 %

RESULTS		
PARAMETER	INITIAL	FINAL
TRANS	1,9 %	1,57 %
DLG	1,7 mm	1,4 mm

**Literature:**

- DLG: 1.9-2.6 mm (Van Esch *et al*,02)
- TRANS: 1.5-1.8 % (Van Esch *et al*,02)

## 2.6. TPS TESTS (Van Esch *et al*, 02)

### QA RAPIDARC TEST PLANS


#### ACCURACY OF DOSE CALCULATION FOR DIFFERENT VOLUMES

**Setup:** 2D array with Octavius

- 2 adjacent target volumes.
- Plan 1 – targets with the same dose
- Plan 2 – targets with different doses

**Analysis:** Gamma analysis (3%, 3mm)

RESULTS		
PLAN	% $\gamma < 1$	Mean $\gamma$
1	99.7 %	0.40
2	99.7 %	0.37



## 2.6. TPS TESTS (Van Esch *et al*, 02)

### QA RAPIDARC TEST PLANS


#### ACCURACY OF DOSE CALCULATION WITH HETEROGENEITIES

**Setup:** Ionization chamber with CIRS phantom (air insert)

- Single target with homogeneous dose.
- Target Sizes: 5x5 cm; 5x10 cm; 5x20 cm.
- Arc: 90° to 270° to exclude couch.

**Analysis:** 3% Tolerance

RESULTS			
PLAN	5x5	5x10	5x20
Dose Difference	2.2 %	2.1 %	2.2 %




## 2.6. TPS TESTS

### CLINICAL RAPIDARC TEST PLANS

**Setup:** 2D array, ionization chamber and EBT2 film with Octavius

- Single arcs exported from hypothetical patient plans onto phantoms.
- Ion chamber (low dose gradients) and 2D array used for **absolute dosimetry**.
- Film used for **relative dosimetry**.



## 2.6. TPS TESTS

### CLINICAL RAPIDARC TEST PLANS

**Analysis:** Gamma analysis (3%, 3mm, 20% threshold)

3% for absolute dosimetry with ion chamber

RESULTS					
PLAN	I.C.	FILM		2D ARRAY	
		% $\gamma < 1$	Mean $\gamma$	% $\gamma < 1$	Mean $\gamma$
1	-0.9 %	97.6 %	0.41	95.4 %	0.41
2	-0.1 %	99.3 %	0.35	97.4 %	0.38
3	1.8 %	99.7 %	0.34	100.0 %	0.34
4	0.1 %	99.4 %	0.33	96.6 %	0.41
5	-1.0 %	98.2 %	0.42	96.6 %	0.45
6	1.7 %	99.9 %	0.38	99.3 %	0.34
<b>AVERAGE</b>	<b>0.3 %</b>	<b>99.0 %</b>	<b>0.37</b>	<b>97.6 %</b>	<b>0.39</b>

**Literature:**  $\leq 3\%$   $\geq 95\%/90\%$   $\geq 95\%$

# 3. Concluding Remarks

---

## MAIN QA TESTS

### ***Daily:***

- Continuous Stripes
- MLC Speed
- RA QA Combined tests

### ***Monthly:***

- Sweeping gap

### ***Patient:***

- Ion chamber
- 2D array
- Film

## MAIN QA CHALLENGES

### ***MLC:***

- Time consuming when done quantitatively to every leaf
- More studies needed with dynalog file viewer analysis

### ***Gantry Speed:***

- No standard tests to evaluate gantry speed (apart from Ling *et al.*,08)

# 3. Concluding Remarks

---

VMAT was implemented successfully in clinical practice.



*Thank you*